

# POREX® High Release Media Flat Sheet

## 36777

Parameters	Value
Technology	Bonded Fiber
Material	HRM Fiber/Polyolefin
Property	Hydrophilic
Shape	Sheet
Width (in)	7.807
Length	N/A
Thickness (in)	0.020
Density (g/cc)	0.250

**Description:**

A 0.5mm T hydrophilic flat sheet with high density HRM (Hydrophilic PE + PP).

Final tolerances for commercial production will be supplied and agreed upon with a *Mutually Understood and Acceptable Specification Document*.

**Disclaimer**

This standard sample is for evaluation of form, fit, and function to determine product specification requirements for custom products. Please reference the values in our data sheet for specification guidelines of standard products. In Vitro Diagnostic (IVD), medical-surgical and pharmaceutical samples are for investigational use only and subject to specific application, test condition and regulatory requirements.

This information is confidential between Porex Corporation (Filtration Group) and the indicated customer.



POREX® Corporation  
500 Bohannon Road  
Fairburn, GA 30213 USA  
info.porex@filtrationgroup.com